


# ILLUSTRATION DRAWING REFERENCE ONLY

REVISIONS			
REV.	DESCRIPTION	DATE	APPROVED
A	INTERNAL CHANGE	5/16/05	SSL
B	5.68±.12 WAS 5.68±.18; ECN 11062	2/23/11	LES



## CUSTOMER APPROVAL

APPROVED BY: \_\_\_\_\_

DATE: \_\_\_\_\_

THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF GAUTHIER BIOMEDICAL, INC., ANY REPRODUCTION IN PART OR WHOLE WITHOUT THE WRITTEN PERMISSION OF GAUTHIER BIOMEDICAL, INC. IS PROHIBITED.

### INSPECTION SYMBOLS



UNLESS OTHERWISE SPECIFIED  
DIMENSIONS ARE IN INCHES  
TOLERANCES ARE:  
FRACTIONS DECIMALS ANGLES  
±1/32 .X ±.020 ±1°  
.XX ±.010  
.XXX ±.005  
.XXXX ±.0005

MATERIAL SEE DETAILS  
HEAT TREAT  
FINISH  
PASSIVATE

**GAUTHIER**   
BIOMEDICAL

GAUTHIER BIOMEDICAL, INC.  
1235 NORTH DAKOTA DRIVE, SUITE G  
GRAFTON, WI 53024  
(262) 546-0010 ph. (262) 546-0011 fax  
E-mail: engineering@gauthierbiomedical.com

APPROVALS	DATE
DRAWN TJM	11/11/04
CHECKED <i>[Signature]</i>	3/03/11
APPROVED <i>[Signature]</i>	3/03/11

DO NOT SCALE DRAWING

DESCRIPTION
MINI 3-SIDED FIXED HANDLE W\SMALL A-O ADAPTOR

SIZE	PART NUMBER	REV.
A	G131150 (cc)	B

SCALE 2:1 CAD FILE: G131150 (cc) SHEET 1 OF 1