


# ILLUSTRATION DRAWING -- REFERENCE ONLY

REVISIONS			
REV.	DESCRIPTION	DATE	APPROVED



THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF GAUTHIER BIOMEDICAL, INC. ANY REPRODUCTION IN PART OR WHOLE WITHOUT THE WRITTEN PERMISSION OF GAUTHIER BIOMEDICAL, INC IS PROHIBITED.

**CUSTOMER APPROVAL**  
 APPROVED BY: \_\_\_\_\_  
 DATE: \_\_\_\_\_

**INSPECTION SYMBOLS**  
 CRITICAL  
 AQL

UNLESS OTHERWISE SPECIFIED  
 DIMENSIONS ARE IN INCHES  
 TOLERANCES ARE:  
 FRACTIONS DECIMALS ANGLES  
 +1/32 .005 .01°  
 .XX ±.010  
 .XXX ±.005  
 .XXXX ±.0005

**gb**  
**Gauthier Biomedical**  
 APPROVALS  
 DRAWN RFM 5/22/09  
 CHECKED RCB 4/2/09  
 APPROVED RCB 6/2/09  
 DO NOT SCALE DRAWING  
 ANODIZE

GAUTHIER BIOMEDICAL, INC.  
 1235 North Dakota Drive, Suite G  
 Grafton, WI 53024  
 (262) 546-0010 ph. (262) 546-0011 fax  
 E-mail: engineering@gauthierbiomed.com

DESCRIPTION  
 OFFSET T-HANDLE TORQUE WR.  
 1/4" SQ. DRIVE

SIZE PART NUMBER  
**B** G403150\_HT (cc) REV. NC  
 SCALE 1:1 CAD FILE: G403150-HT (cc) SHEET 1 OF 1