


ILLUSTRATION DRAWING REFERENCE ONLY

REVISIONS			
REV.	DESCRIPTION	DATE	APPROVED
A	ECN 8177	6/26/08	GDD



CUSTOMER APPROVAL

APPROVED BY: _____
DATE: _____

THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF GAUTHIER BIOMEDICAL, INC., ANY REPRODUCTION IN PART OR WHOLE WITHOUT THE WRITTEN PERMISSION OF GAUTHIER BIOMEDICAL, INC. IS PROHIBITED.

INSPECTION SYMBOLS

-  CRITICAL
-  AQL

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES
TOLERANCES ARE:
FRACTIONS DECIMALS ANGLES
 $\pm 1/32$.X $\pm .020$ $\pm 1^\circ$
.XX $\pm .010$
.XXX $\pm .005$
.XXXX $\pm .0005$



Gauthier Biomedical, Inc.

GAUTHIER BIOMEDICAL, INC.
1235 NORTH DAKOTA DRIVE, SUITE G
GRAFTON, WI 53024
(262) 546-0010 ph. (262) 546-0011 fax
E-mail: engineering@gauthierbiomedical.com

MATERIAL SEE DETAILS

HEAT TREAT

FINISH

PASSIVATE

APPROVALS	DATE
DRAWN MCS	1/23/08
CHECKED <i>MCS</i>	6/26/08
APPROVED <i>MCS</i>	6/26/08
DO NOT SCALE DRAWING	
ANODIZE	

DESCRIPTION
TORQUE LIMITER
BALANCED SILICONE T-HANDLE
W/BLANK EXTRA LONG SHAFT

SIZE PART NUMBER
A G412154 (CC) REV. A
SCALE 0.75:1 CAD FILE: G412154 (cc) SHEET 1 OF 2


ILLUSTRATION DRAWING REFERENCE ONLY

REVISIONS			
REV.	DESCRIPTION	DATE	APPROVED
A	ECN 8177	6/26/08	GDD

17-4 STAINLESS
TOP CAP



17-4 STAINLESS
BOTTOM CAP



300 SERIES STAINLESS
END CAPS

CUSTOMER APPROVAL

APPROVED BY: _____

DATE: _____

THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF GAUTHIER BIOMEDICAL, INC.. ANY REPRODUCTION IN PART OR WHOLE WITHOUT THE WRITTEN PERMISSION OF GAUTHIER BIOMEDICAL, INC. IS PROHIBITED.

INSPECTION SYMBOLS

 CRITICAL

 AQL

UNLESS OTHERWISE SPECIFIED
DIMENSIONS ARE IN INCHES
TOLERANCES ARE:
FRACTIONS DECIMALS ANGLES
±1/32 .X ±.020 ±1°
.XX ±.010
.XXX ±.005
.XXXX ±.0005



Gauthier Biomedical, Inc.

GAUTHIER BIOMEDICAL, INC.
1235 NORTH DAKOTA DRIVE, SUITE G
GRAFTON, WI 53024
(262) 546-0010 ph. (262) 546-0011 fax
E-mail: engineering@gauthierbiomedical.com

MATERIAL SEE DETAILS

HEAT TREAT

FINISH

PASSIVATE

APPROVALS

DRAWN

CHECKED

APPROVED

DO NOT SCALE DRAWING

ANODIZE

DATE

1/23/08

6/26/08

6/26/08

DESCRIPTION

TORQUE LIMITER
BALANCED SILICONE T-HANDLE
W/BLANK EXTRA LONG SHAFT

SIZE

A

PART NUMBER

G412154 (CC)

REV.

A

SCALE 1:1

CAD FILE: G412154 (CC)

SHEET 2 OF 2